Pain-Free Melting Point Determination

James Lee, Ph.D. Stanford Research Systems



Who is Stanford Research Systems?



- In business since 1980
- Full catalog is over 200 pages
- Famous for first digital lock-in amplifier
- Successful quadrupole mass spectrometer (RGA)
- Now makes 3 kinds of melting point apparati



135.4 °C



Introduction



Student grade melting point apparatus Integral RTD thermometer Microprocessor controlled temperature ramps PID gives fast preheats without overshoot Tube Tapper (integral capillary packing device) Easy to clean, maintain <u>Price (single unit)</u> : \$790

What is an RTD?

Resistance Temperature Detectors (RTDs) are simply materials whose resistivity changes as a function of temperature.



What is PID?

P.I.D. control allows tight temperature control of the oven using 3 types of error correction:

Proprotional – Make a correction term proportional to the error

Integral – Make a correction term that integrates the error with respect to time

Derivative – Make a correction term that differentiates the error with respect to time

Correction = $Pe(t) + I \int e(t) dt + D de/dt$ (where e(t) is the error signal with respect to time, P, I, and D are constants)

Temperature Measurement

MEL-TEMP

Liquid in-glass (often Mercury) Accuracy typically +/- 2 degree Resolution typically 0.2 degree

DigiMelt

RTD accuracy is typically ± 0.6 degree C <200 C ± 1.0 degree C >200 C

Resolution is 0.1 degree C



Note that the accuracy of the thermometer isn't really the accuracy of the MelTemp. The device must be calibrated against melting point standards to give a slope and offset correction to the thermometer reading.

SRS takes care of calibration for you. We calibrate every unit we sell to our stated accuracy.



Temperature Measurement

MEL-TEMP

500 C is reachable, thermometer often only goes to 400 C

Thermometer can easily be broken or taken!

DigiMelt

Temperatures > 260 C disallowed

RTD is never missing

RTD is difficult to destroy





Temperature Control

MEL-TEMP

Open Loop Variac: No feedback Student is the controller Typical student is a poor controller Overshoot leads to fanning the block which leads to breakage



DigiMelt

Full PID control with RTD feedback

Microprocessor lets students focus on their samples

Ramp rates of 0.5, 1, 2, 5 C

Oven turns off after 30 minutes of idle



MelTemp Delivers a Set Power, not a Set Temperature



When students turn up the knob on the variac, they often don't realize that the delivered ramp is very non-linear with the power setting. This often leads to students turning up the variac too high and overshooting their target temperature. This wastes time and frustrates many students.

DigiMelt Follows a User-Defined Temperature Program



Time

Sample Viewing

MEL-TEMP

Hard to see all three tubes

Light reflects from sample to eye via window, lens

Student must switch between viewing sample and viewing thermometer



DigiMelt

Large lens allows view of three samples simultaneously

Student records critical temperatures (onset, meniscus, clear point) by touching keypad





Using DigiMelt : Enter Settings



The keypad makes setting up a melt program very easy. The students push the yellow 'START TEMP' and then use the blue buttons with arrow keys to set the starting temperature. The students do the same for the 'RAMP RATE' and 'STOP TEMP' keys. Once the parameters are entered, they push the 'START/STOP' button and the microprocessor takes care of heating the block.

Using DigiMelt to Record Data









Students can 'flag' critical points in the melting process with the blue buttons on the keypad. Above are views through the lens showing a sample of vanillin during pre-heat, at the meniscus point, and at the clear point. After the melt is over (cooling LED is lit), the data are retrived by pushing the buttons again.

Packing Capillaries



When students tap sample into capillaries, lots of capillaries end up on the floor, lots of capillaries are broken, wasted

DigiMelt

Integral cell phone vibrator motor "Tube Tapper" Capillaries can be packed 3 at a time Students can get 3 samples with the same sample height



MEL-TEMP

Incandescent bulb will burn out

Mean time between failures is about 1000 hrs

Broken thermometer = \$200 (25%)

DigiMelt

White LED has >100,000 hours mean time between failures

RTD is well-protected





MEL-TEMP

Capillary holder uses socket head cap screws or other screws



DigiMelt No tools required



MEL-TEMP

Capillary holder uses socket head cap screws or other screws



DigiMelt No tools required



MEL-TEMP

Lens, window held with snap rings



DigiMelt No tools required





Safety Certification



Student safety is paramount in our design. SRS had the entire DigiMelt design independently verified to comply with UL safety regulations by MET labs. Every unit undergoes safety testing to maintain our certification.

Summary

MEL-TEMP

Students know there's a better way

Wasted time Hard to use equipment Broken thermometers

FRUSTRATION

DigiMelt

Students focus on the sample: THE REAL EXPERIMENT Faster labs Safer labs Less waste Lower cost of ownership



