



**Stanford Research Systems, Inc.**  
290 C Reamwood Avenue, Sunnyvale, CA 94089, U.S.A.  
Phone: 408-744-9040 FAX: 408-744-9049  
Email: [janied@thinkSRS.com](mailto:janied@thinkSRS.com) .. Website: [www.thinkSRS.com](http://www.thinkSRS.com)

---

10/8/2019

### ***Volatility Statement for SR785***

The SR785 contains various types of memories, as listed below:

#### **Non-Volatile Memory:**

(2KB x 8) x 3 EPROM      Not user modifiable;  
Contains instrument's firmware and factory  
calibration bytes

#### **Volatile Memory:**

(128KB x 5) x 2 SRAM      user modifiable; used for user settings.  
To clear:  
Depress the backspace key while powering ON the unit

(128KB x 6) x 2 SRAM      user modifiable; used for data storage  
Depress the backspace key while powering ON the unit

Please feel free to contact SRS should you have any further questions.

Janie Du  
Sales & Marketing