



DETAIL A
SCALE 8.000
DIMENSIONS ARE TYPICAL



REVISIONS			
REV	DESCRIPTION	DATE	APPROVED
B	CHANGED INSIDE HOLE DIAMETER CHANGED TOLERANCE.	4/25/22	GAW



DETAIL B
SCALE 4.000
3 EQUAL SPACES

SEE DETAIL B

NOTES: UNLESS OTHERWISE SPECIFIED

- 1 MATERIAL: STAINLES STEEL COIL 302/304 HALF HARD.
.010 THICK.
- 2 FINISH: BLACK OXIDE.
- 3. REMOVE ALL BURRS AND SHARP EDGES.
- 4. PART MARKING: BAG AND TAG WITH VENDOR ID,
MONTH/YEAR OF MANUFACTURE PART NUMBER AND
CURRENT REVISION LEVEL.

MATERIAL	UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES.	STANFORD RESEARCH SYSTEMS, INC. 1290-D Redwood Avenue, Sunnyvale CA 94089
FINISH		
TOLERANCES:	.x ±	DRAWING TITLE: CHOPPER BLADE, 25-30 SLOT ENCODER DISC, SR542
	.xx ±	
	.xxx ± .005	
ANGLE ±		
ENG BY/DATE A BERGER	SIZE/PROVE FILENAME: D	SRS PART NO. 7-02828-716
DWN BY/DATE G.WATERS	Jun-04-20	SCALE: 4.00
DO NOT SCALE DRAWING		SHEET: 1 OF 1

NEXT ASSY	USED ON
-----------	---------