

| Certificate of Volatility | | | | |
|---|------------------------|--|--|---|
| Model: PS310, 325, 350 | Part Number: | Manufacturer: Stanford Research Systems, Inc | | |
| | | Street Address: 1290-D Reamwood Ave | | |
| | | City: Sunnyvale | State: CA | Zip: 94089 |
| Item Nomenclature: | | | | |
| Volatile Memory | | | | |
| Does the item contain volatile memory (i.e., memory whose contents are lost when power is removed)? <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | | | | |
| If the answer is 'Yes', please provide the following information for each type (use additional sheets if required): | | | | |
| Type (SRAM, DRAM, etc.): | Size: | User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No | Function: | Process to Clear: |
| Type (SRAM, DRAM, etc.): | Size: | User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No | Function: | Process to Clear: |
| Type (SRAM, DRAM, etc.): | Size: | User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No | Function: | Process to Clear: |
| Non-Volatile Memory | | | | |
| Does the item contain non-volatile memory (i.e., memory whose contents are retained when power is removed)? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | | | | |
| If the answer is 'Yes', please provide the following information for each type (use additional sheets if required): | | | | |
| Type (BBRAM, Flash, EEPROM, etc.): EEPROM | Size: 32 kB | User Modifiable: <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | Function: Firmware storage and factory calibration | Process to Clear: Programmed at factory only |
| Type (BBRAM, Flash, EEPROM, etc.): Battery backed up SRAM | Size: 2 kB | User Modifiable: <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No | Function: Stores working calibration data and stored configurations | Process to Clear: Holding down the CLR button while turning on power will erase all stored configurations and reload the factory calibration data. Or remove the back-up battery. |
| Type (BBRAM, Flash, EEPROM, etc.): | Size: | User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No | Function: | Process to Clear: |
| Media | | | | |
| Does the item contain media storage capability (i.e., removable or nonremovable disk drives, tape drives, memory cards, etc.)? <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No | | | | |
| If the answer is 'Yes', please provide the following information for each type (use additional sheets if required): | | | | |
| Type (Disk, Tape, etc.): | Size: | User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No | Function: | Process to Clear: |
| Removable: <input type="checkbox"/> Yes <input type="checkbox"/> No | | | | |
| Type (Disk, Tape, etc.): | Size: | User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No | Function: | Process to Clear: |
| Removable: <input type="checkbox"/> Yes <input type="checkbox"/> No | | | | |
| Type (Disk, Tape, etc.): | Size: | User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No | Function: | Process to Clear: |
| Removable: <input type="checkbox"/> Yes <input type="checkbox"/> No | | | | |
| Additional Information: | | | | |
| | | | | |
| Vendor/Manufacturer Representative Information | | | | |
| Name: JANIE DU | Title: SALES | Office Phone: 408-744-9040 | Fax/Email: INFO@THINKSRS.COM | |