

Certificate of Volatility				
Model: SIM984	Part Number:	Manufacturer: Stanford Research Systems		
		Street Address: 1290-D Reamwood Ave		
		City: Sunnyvale	State: CA	Zip: 94089
Item Nomenclature:				
Volatile Memory				
Does the item contain volatile memory (i.e., memory whose contents are lost when power is removed)? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No				
If the answer is 'Yes', please provide the following information for each type (use additional sheets if required):				
Type (SRAM, DRAM, etc.): SRAM	Size: 1kbyte	User Modifiable: <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	Function: program varibal storage and communications data buffer	Process to Clear: turn off power
Type (SRAM, DRAM, etc.):	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Type (SRAM, DRAM, etc.):	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Non-Volatile Memory				
Does the item contain non-volatile memory (i.e., memory whose contents are retained when power is removed)? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No				
If the answer is 'Yes', please provide the following information for each type (use additional sheets if required):				
Type (BBRAM, Flash, EEPROM, etc.): Flash	Size: 32kbyte	User Modifiable: <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	Function: executable progam storage	Process to Clear: desolder microcontroller -- destructive to product
Type (BBRAM, Flash, EEPROM, etc.): EEPROM	Size: 768 bytes	User Modifiable: <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	Function: stores factory info & last user setting	Process to Clear: Set instrument to defaults: Gain 1, Bandwidth 100 Hz. Then remove power
Type (BBRAM, Flash, EEPROM, etc.):	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Media				
Does the item contain media storage capability (i.e., removable or nonremovable disk drives, tape drives, memory cards, etc.)? <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No				
If the answer is 'Yes', please provide the following information for each type (use additional sheets if required):				
Type (Disk, Tape, etc.): Removable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Type (Disk, Tape, etc.): Removable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Type (Disk, Tape, etc.): Removable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Size:	User Modifiable: <input type="checkbox"/> Yes <input type="checkbox"/> No	Function:	Process to Clear:
Additional Information:				
Vendor/Manufacturer Representative Information				
Name: Janie Du	Title: Sales	Office Phone: 408-744-9040	Fax/Email: 408-744-9049/ info!thinkSRS.com	