

## ***Volatility Statement for the SRS SIM964***

The SIM964 contains various types of memories, as listed below:

### **Non-Volatile Memory:**

|                    |  |
|--------------------|--|
| 32 KB Flash EEPROM | Not user modifiable; contains instrument's firmware  |
| 768 Bytes EEPROM   | user modifiable; contains calbytes & instrument settings<br><u>To clear instrument settings:</u><br>Send '*RST' and power cycle the unit.<br>Note: calbytes are not user-accessible. |

### **Volatile Memory:**

|          |   |
|----------|---|
| 1 KB RAM | user modifiable; used for operation and user settings.<br><u>To clear:</u><br>Power cycle the unit. |
|----------|---|

Please feel free to contact SRS should you have any further questions.

Stanford Research Systems  
1290-D Reamwood Avenue  
Sunnyvale, CA 94089

(408)744-9040  
[www.thinksrs.com](http://www.thinksrs.com)