




REVISIONS			
REV	DESCRIPTION	DATE	APPROVED
B	CHANGED TWO CNSK HOLES TO .230 DIA	08/01/88	AWR




TOP VIEW



FRONT VIEW




RIGHT SIDE VIEW



- NOTES: 1) ALL DIMENSIONS ARE IN INCHES
2) TOL +/- .010
3) USE MINIMUM RADII ON ALL BENDS
4) MATERIAL: .090 ALUMINUM
5) FINISH: BLACK ANODIZE
DO NOT GRAIN, TUMBLE O.K.

STANFORD RESEARCH SYSTEMS, Inc.					
1290-D REAMWOOD AVENUE SUNNYVALE, CA 94089 PHONE: (408) 744-9040 FAX: 4087449049 TELEX: 706891 SRS UD			ASSEMBLY: SR510/PS300 SCIENTIFIC INSTRUMENTS		
			TITLE: RACK MOUNT (TWO PER ASSEMBLY)		
DRAWN BY: AWR	SCALE: 1:1	DATE: 03/11/88	SHEET 1 OF 1	SIZE C	REV B
				SRS REF# SR510-26	SRS PART# 7-00205-720


REVISIONS			
REV	DESCRIPTION	DATE	APPROVED
C	ADDED TWO (2) #10-32 TAPPED HOLES	12/19/94	JG



- NOTES: UNLESS OTHERWISE SPECIFIED
- 1) ALL DIMENSIONS ARE IN INCHES AND APPLY AFTER PLATING
 - 2) TOLERANCE: +/- .005 UNLESS OTHERWISE NOTED
 - 3) MATERIAL: 6061 ALUMINUM
 - 4) DOVETAILS TO BE MACHINED FOR SMOOTH, SLIDING FIT OF TWO (2) MATING PARTS
 - 5) FINISH: ANODIZE, BLACK

STANFORD RESEARCH SYSTEMS, Inc.			
1290-D REAMWOOD AVENUE SUNNYVALE, CA 94089		ASSEMBLY: HIGH VOLTAGE POWER SUPPLY	
PHONE: (408) 744-9040 FAX: 4087449049 TELEX: 706891 SRS UD		TITLE: RACKMOUNT, DOUBLE	
SIZE C	SRS REF# SR300-38	SRS PART# 7-00194-720	REV C
DRAWN BY: LD	SCALE: 2:1	DATE: 7-29-88	SHEET 1 OF 1

REVISIONS			
REV	DESCRIPTION	DATE	APPROVED



- NOTES: 1) ALL DIMENSIONS ARE IN INCHES
2) TOL +/- .010
3) REMOVE ALL SHARP EDGES AND BURRS (GRAINING O.K.)
4) MATERIAL: .125" THICK (8 ga.) ALUMINUM STRIP
5) FINISH: BLACK ANODIZE
6) RUBBER OR METAL STAMP ON THIS PART (LOCATION OPTIONAL) IN .06 - .12 HIGH CHARACTERS: 7-00676-720.REVO

STANFORD RESEARCH SYSTEMS, Inc.			
1290-D REAMWOOD AVENUE SUNNYVALE, CA 94089 PHONE: (408) 744-9040 FAX: 4087449049 TELEX: 706891 SRS UD		ASSEMBLY: DOUBLE RACK MOUNT - PS,DS34X,SR560/70,SR630	
		TITLE: SAFETY LOCKING PLATE	
DRAWN BY: JG		SIZE: C SRS REF# PS300-52	SRS PART# 7-00680-720 REV 0
SCALE: 2:1		DATE: 12/19/94	SHEET 1 OF 1

Use this kit to mount two half-rack sized PS310, PS325, or PS350 instruments side-by-side into a standard 19-inch rack.
Included hardware: 2x Side mounting brackets, 2x Rack handles, 1x Center connecting latch, 1x Safety locking plate, and screws.

