




REVISIONS			
REV	DESCRIPTION	DATE	APPROVED
B	CHANGED TWO CNSK HOLES TO .230 DIA	08/01/88	AWR




TOP VIEW



FRONT VIEW



RIGHT SIDE VIEW



- NOTES: 1) ALL DIMENSIONS ARE IN INCHES
2) TOL +/- .010
3) USE MINIMUM RADII ON ALL BENDS
4) MATERIAL: .090 ALUMINUM
5) FINISH: BLACK ANODIZE
DO NOT GRAIN, TUMBLE O.K.

STANFORD RESEARCH SYSTEMS, Inc.					
1290-D REAMWOOD AVENUE SUNNYVALE, CA 94089 PHONE: (408) 744-9040 FAX: 4087449049 TELEX: 706891 SRS UD			ASSEMBLY: SR510/PS300 SCIENTIFIC INSTRUMENTS		
			TITLE: RACK MOUNT (TWO PER ASSEMBLY)		
			SIZE C	SRS REF# SR510-26	SRS PART# 7-00205-720
DRAWN BY: AWR			SCALE: 1:1	DATE: 03/11/88	SHEET 1 OF 1



Use this kit to mount one half-rack sized SRS instrument into a standard 19-inch rack.

Included hardware: 1x Side mounting bracket, 1x Extension panel, 2x Rack handles, and screws.