



DETAIL A
SCALE 8.000
DIMENSIONS ARE TYPICAL

- NOTES: UNLESS OTHERWISE SPECIFIED
- 1 MATERIAL: STAINLES STEEL COIL 302/304 HALF HARD. .010 THICK.
 - 2 FINISH: BLACK OXIDE.
 3. REMOVE ALL BURRS AND SHARP EDGES.
 4. PART MARKING: BAG AND TAG WITH VENDOR ID, MONTH/YEAR OF MANUFACTURE PART NUMBER AND CURRENT REVISION LEVEL.
 5. PART TO BE RoHS COMPLIANT.



DETAIL B
SCALE 4.000
3 EQUAL SPACES

SEE DETAIL B

REVISIONS			
REV	DESCRIPTION	DATE	APPROVED
B	OD FOR 10 SLOT TRACK CHANGED TO Ø3.050.	5/28/20	AB
C	CHANGED ID AND TOLERNACE ADDED RoHS COMPLIANT NOTE.	2/14/22	GAW
D	CHANGED MATERIAL SPECIFICATION	3/22/22	GAW

□ .010

MATERIAL	UNLESS OTHERWISE SPECIFIED, DIMENSIONS ARE IN INCHES.	STANFORD RESEARCH SYSTEMS, INC. 1290-D Reamwood Avenue, Sunnyvale CA 94089	
FINISH			
1	TOLERANCES:	DRAWING TITLE:	
2	.x ±	CHOPPER BLADE, 10-100 SLOT ENCODER DISC, SR542	
	.xx ±	SIZE	PROVE FILENAME:
	.xxx ± .005	D	
	ANGLE ±	SRS PART NO.	REV
ENG BY/DATE	A BERGER	7-02741-716	D
DWN BY/DATE	G.WATERS Jun-04-20	SCALE: 4.00	SHEET: 1 OF 1
NEXT ASSY	USED ON	DO NOT SCALE DRAWING	