



**Stanford Research Systems, Inc.**  
290 C Reamwood Avenue, Sunnyvale, CA 94089, U.S.A.  
Phone: 408-744-9040 FAX: 408-744-9049  
Email: [janied@thinkSRS.com](mailto:janied@thinkSRS.com) .. Website: [www.thinkSRS.com](http://www.thinkSRS.com)

---

## **Volatility Statement for SR810 and SR830**

**SRAM, 256 k bytes** is used for settings, variables, and data buffers.

This SRAM is supplied by a 3.3V battery when the unit is not connected to the power mains.

Holding the **[Setup]** key during power-on will restore the settings and variables to their factory default values.

**Sanitizing the data buffers requires over-writing them with zero (or Any non-sensitive) data:**  
(The following assumes the unit has already been set to factory defaults) Disconnect the input.

Set the remote interface to match your computer interface (RS232 or GPIB).

Query \*IDN? to verify that the remote connection is working.

Send the following two commands. The first sets the data recording rate to maximum (500 Hz). The second starts recording.

```
SRAT 13  
STRT
```

Wait 35 seconds. The unit will make a "done" sound sequence. It is now completely cleaned of stored data.

Please feel free to contact SRS should you have any further questions.